Government of India Atomic Energy Regulatory Board Radiological Safety Division Niyamak Bhavan, Anushaktinagar, Mumbai-400094

FORMAT FOR QUALITY ASSURANCE TEST FOR DIAGNOSTIC X-RAY EQUIPMENT

(Applicable for Radiography/Radiography&Fluoroscopy/C-Arm/InterventionalRadiology Equipment)

(Periodic Quality Assurance shall be carried out at least once in two years and also after any repairs having radiation safety implications)

A. DETAILS OF THE DIAGNOSTIC X-RAY EQUIPMENT

1	Name of the Institution and City	
2	Type of Equipment	
3	Model Name	
4	Name of the Manufacturer	
5	Name(s) of Person(s) testing the equipment	
	and Name of Supplier/Service Agency	
6	Dates and Duration of the Tests	

B. SUMMARY OF RADIATION SAFETY PERFORMANCE TEST REPORT:

Sr.	Parameters	Specified	Measured	Teleronce	Remarks
No.	Tested	Values	Values	Tolerance	
1	Congruence of radiation and optical field	Measurement at 60 kV		Tolerance : I X I + I X' I \leq 2% of FFD I Y I + I Y' I \leq 2% of FFD	
2	Central Beam Alignment	Measurement at 60 kV		Central Beam Alignment $< 1.5^{\circ}$	
3	Effective focal spot measurement FFD= 60 cm			Tolerance : 1. + 0.5 f for f < 0.8 mm 2. + 0.4 f for $0.8 \le f \le 1.5$ mm 3. + 0.3 f for f > 1.5 mm	
4	Accuracy of Operating Potential (kV)	(last value maximum kV)		± 5 kV	
5	Accuracy of Irradiation Time(sec.)			% Error < 10 %	
6	Total filtration	Measurement at maximum kV		Tolerance : 1.5 mm Al for $kV \le 70$ 2.0 mm Al for $70 < kV \le 100$ 2.5 mm Al for $kV > 100$	
7	Linearity of mA/mAs loading Stations			CoL<0.1	

8	Consistency		$CoV \leq 0.05$	
	of radiation			
	output			
9	Low contrast	Measurement	3.0 mm hole pattern must be	
	resolution	shall be done	resolved	
		on each mode		
		of Image		
		Intensifier		
10	High contrast	Measurement	1.5 lp/mm pattern must be	
	resolution	shall be done	resolved	
		on each mode		
		of Image		
		Intensifier		
11	Exposure Rate	Applicable	Tolerance :	
	at Tabletop	Fluoroscopy	1. Exposure Rate without AEC	
		Mode	$mode \le 5 cGy/Min$	
			2. Exposure Rate with AEC	
			mode $\leq 10 \text{ cGy/Min}$	
12	Radiation	Measurement	Tube Leakage < 1 mGy in one	
	leakage level	at Maximum	hour	
	at 1m from	kVp and		
	tube housing	corresponding		
	and	tube current		
	Collimator			
	(Pl. mention			
	values			
	separately)			

Note:

- 1. Tests as mentioned in Sr. No. 1,2,3 are not applicable for C-Arm/Interventional Radiology Equipment.
- 2. Tests as mentioned in Sr. No. 9,10,11 are not applicable for only radiography mode of x-ray equipment.

I hereby undertake that all the information provided above is correct and in accordance with the detailed Quality Assurance Report enclosed herewith.

Place: Date: Signature: Name of the Service Engineer: Name of Supplier/Service Agency: Seal of Supplier/Service Agency:

#Signature of Institution's Representative: Name of Institution: Seal of the Institution:

Quality Assurance Tests Report shall be signed by Institution's Representative and duly stamped by the User's Institution.

1. CONGRUENCE OF RADIATION & OPTICAL FIELD

Operating parameters:

FFD (cm)	100	kV	60	mAs	

Shift in the edges of the	radiation field			
Dimensions (cm)	Observed shift	% of FFD	Tolerance	Remark
I X I + I X' I			2 % of FFD	
I Y I + I Y' I			2 % of FFD	

2.CENTRAL BEAM ALIGNMENT

Operating parameters:

$\Gamma D (CIII) = 100 KV = 00 IIIAS$

Observe the images of the two steel balls on the radiograph and evaluate tilt in the central beam.Observed tiltRemarkTolerance: Central Beam Alignment< 1.5°</td>

3.EFFECTIVE FOCAL SPOT MEASUREMENT

FFD (cm) 60

	Stated value of Focal	Measured value of Focal Spot	Tolerance:
	Spot	Size (mm x mm)	1. $+ 0.5$ f for f < 0.8 mm
	Size (mm x mm)		2. + 0.4 f for $0.8 \le f \le 1.5$ mm
Large			3. $+ 0.3$ f for f > 1.5 mm
Focus			
Small			
Focus			

4. ACCURACY OF OPERATING POTENTIAL

- **5. ACCURACY OF IRRADIATION TIME**
- 6. TOTAL FILTRATION

FFD (cm) 100

Applied	Set		Measured values							
kVp	time		mA stations							
	(s)	Minim	um mA	mA Routinely used Maximum mA					Average	Remarks
		station		mA	station	ion station			Time (s)	
		kVp	Time (s)	kVp	Time (s)	kVp	Time (s)			
Max.	Max.									
kVp	kVp									
Tolerance for $kVp: \pm 5 kV$										
Tolerance	e for Tim	er: %	Error: ± 10	%						
Total Filt	ration (m	easurem	ent at maxii	num kVp)):	n	nm of Al			
Tolerance	e for Tota	al Filtrat	ion: 1.5 mm	Al for k	$V \le 70, 2.0 \text{ n}$	nm Al fo	$r kV \le 100$,	2.5 mm Al	for $kV > 10$	0

7. LINEARITY OF (mA/mAs) LOADING STATIONS

Operating parameters:

FFD (cm) 100 kV Time (s)		1.0.0		· _ · _ · _ · _ · _ · _ · _ · _ ·	
	FFD (cm)	100	kV	Time (s)	

mA applied	Radi	ation Output (μGy)	Average Output	μGy /mAs (X)	Coefficient of Linearity (CoL)	Remarks
	Reading 1	Reading 2	Reading 3	(µGy)			
Tolerance	e: COL < 0.1						

8. OUTPUT CONSISTANCY

FFD (cm) 100

Applied mAs Radiation Output (µGy)							Average	Coefficient of Remar	Remarks		
ΚV		1	2	3	4	5	(X)	Variation			
								(CoV)			
Tole	ance : (CoV < 0.05									

9. LOW CONTRAST SENSITIVITY:

Diameter of the smallest size hole clearly resolved on the monitor	
Recommended performance standard	3.0 mm hole pattern must be resolved

10. HIGH CONTRAST SENSITIVITY

|--|

Recommended performance standard	1.5 l	p/mm	pattern must be resolved
	1.0 1	P/	

11. EXPOSURE RATE AT TABLE TOP

Distance from Focus	Applied kV	Applied mA	Exposure rate at table top (cGy/Min)	Remarks	
to Table top					
Tolerance :					
1. Exposure Rate without AEC mode \leq 5 cGy/Min					
2. Exposure Rate with AEC mode $\leq 10 \text{ cGy/Min}$					

Minimum focus to table top distance shall be 30 cm for fluoroscopy

12.TUBE HOUSING LEAKAGE

Operating parameters:

FDD (cm)100kVp(Max)mATime(s) \approx 1 sec
--

Location	Exposure level (mR/hr)					Result
(at 1.0 m from the focus)	Left	Right	Front	Back	Тор	
Tube						
Collimator						
Tolerance: Maximum leakage radiation level at 1 meter from the focus should be $\leq 1 \text{ mGy}(114 \text{ mR})$ in one hour						

Work load 180 mAmin in one hr

Max leakage = <u>180 mAmin in 1 hr X ---- Max Exposure level (mR/hr)</u> 60 X -----mA used for measurement

Maximum radiation leakage from tube housing = ------ mR in one hour Maximum radiation leakage from tube collimator = ------ mR in one hour

Details of **Radiation Protection Survey of the installation**

Date of radiation protection survey:

Whether radiation survey meter used for the survey has valid calibration certificate: Yes/No

Equipment Setting:-	
Applied Current (mA):	
Applied Voltage (kV):	
Exposure time(s):	
Workload:	
Provide the measured maximum radiation levels (mR/hr) at different locations	

Control console(Operator Position)	
Outside patient entrance door	
Behind Windows (if applicable)	
Patient Waiting Area	

Maximum Radiation level/week (mR/wk) = <u>---- mAmin/week X ----Max. radiation level (mR/hr)</u> 60 X -----mA used for measurement

Permissible limit

For location of Radiation Worker: 20 mSv in a year (40 mR/week) For Location of Member of Public: 1 mSv in a year (2mR/week)