

PERFORMANCE TEST REPORT FOR DENTAL X-RAY EQUIPMENT

(Periodic Quality Assurance shall be carried out at least once in two years and also after any repairs having radiation safety implications)

A. DETAILS OF THE DIAGNOSTIC X-RAY EQUIPMENT

1	Name of the Institution and City	
2	Type of Equipment	
3	Model Name	
4	Name of the Manufacturer	
5	Name(s) of Person(s) testing the equipment and Name of Supplier/Service Agency	
6	Dates and Duration of the Tests	

B. SUMMARY OF RADIATION SAFETY PERFORMANCE TEST REPORT:

Sr. No.	Parameters Tested	Measured Values/Observation	Tolerance	Remark
3.1	Operating Potential		± 5 kV	
3.2	Accuracy of Irradiation time (second)		% Error $\leq 10\%$	
3.3	Minimum Total Filtration		Min. 1.5 mm Al for kVp < 70 Min. 2.0 mm Al for 70 < kVp < 100 Min. 2.5 mm Al for kVp > 100	
3.4	Timer linearity (Coefficient of Linearity)		< 0.1	
3.5	mA/mAs Linearity		CoL < 0.1	
3.6	Consistency of Radiation Output (Coefficient of Variation)		< 0.05	
3.7	Radiation leakage level from X-ray tube housing at 1m from the focus (measured at max kVp) (mGy in one hour)		Tube Leakage < 0.25 mGy in one hour (For Dental Intraoral) Tube Leakage < 1 mGy in one hour (For Dental Extra-oral)	

C. Additional performance test to be carried out for Dental CBCT

Sr. no.	Parameters	Observation (Write "Yes/No/ or give values of the measurement", as the case may be)	Tolerance Limit	Remark
---------	------------	---	-----------------	--------

4.1	Low Contrast Sensitivity		As per technical specification	
4.2	High Contrast Sensitivity		As per technical specification	
4.3	CT dose index (CTDI _w) (*If manufacture specify CTDI _w value then perform this test)		±20% of quoted value	

I hereby undertake that all the information provided above is correct and in accordance with the detailed Quality Assurance Report enclosed herewith.

Place:
Date:

Signature:
Name of the Service Engineer:
Name of Supplier/Service Agency:
Seal of Supplier/Service Agency:

#Signature of Institution's Representative:

Name of Institution:
Seal of the Institution:

Quality Assurance Tests Report shall be signed by Institution's Representative and duly stamped by the User's Institution

1. ACCURACY OF OPERATING POTENTIAL/ACCURACY OF TIMER

*FDD(cm)	
----------	--

Applied kVp	Set time (s)	Measured values						
		mA stations						
		mA station -1		mA station -2		Ave KVp	Ave Time (s)	Remarks
		kVp	Time (s)	kVp	Time (s)			
Max. kVp								
Tolerance for kVp : ± 5 kV								
Tolerance for irradiation time : ± 10 %								
Measured Minimum Total Filtration :mm of Al at maximum kVp								
Tolerance : 1.5 mm Al for kV ≤ 70, 2.0 mm Al for kV ≤ 100, 2.5 mm Al for kV > 100								

**FDD stands for focal spot to detector distance*

2.LINEARITY OF mA LOADING STATION

Operating parameters

FDD(cm)	100	kV		Time	
---------	-----	----	--	------	--

mA Range	Output (μGy)	Average (μGy)	μGy /mAs (X)	Coefficient of Linearity (CoL)	Remarks

	Reading 1	Reading 2				

Tolerance: CoL < 0.1

3. REPRODUCIBILITY OF RADIATION OUTPUT

FDD(cm)	
---------	--

Applied kV	mAs	Output (µGy)			Average (X)	Coefficient of Variation (CoV)	Remarks
		1	2	3			

Tolerance : COV < 0.05

4. LOW CONTRAST SENSITIVITY:

Diameter of the smallest size hole clearly seen on the monitor	
Recommended performance standard	3 mm hole of measuring test tool must be resolved

5. HIGH CONTRAST SENSITIVITY

Bar strips resolved on the monitor (lp/mm)	
Recommended performance standard	1.5 lp/mm of measuring test tool must be resolved

6. TUBE HOUSING LEAKAGE

Operating parameters:

Distance from the focus (cm)	100	kVp	(Max)	mA		Time(s)	
------------------------------	-----	-----	-------	----	--	---------	--

Location (at 1.0 m from the focus)	Exposure level (mGy/hr or mR/hr)					Result
	Left	Right	Front	Back	Top	
Tube						... mGy in 1 hr
Collimator						... mGy in 1 hr

Tolerance: Maximum leakage radiation level at 1 meter from the focus shall be ≤ 1 mGy in one hour.

Work load = 180 mAmin in one hour (for dental OPG and dental CBCT equipment)

20 mAmin in one hour for dental Intra oral equipment

Max leakage = $\frac{\text{---mAmin in 1 hr} \times \text{---Max leakage (mGy/hr or mR/hr)}}{60 \times \text{----mA used for measurement}}$

Maximum radiation leakage from tube = ----- mGy in one hour

7 Details of Radiation Protection Survey of the installation

Date of radiation protection survey:

Whether radiation survey meter used for the survey has valid calibration certificate: Yes/No

Equipment Setting:-

Applied Current (mA):

Applied Voltage (kV):

Exposure time(s):

Workload:

Provide the measured maximum radiation levels (mR/hr) at different locations

Location	Max. Radiation level (mR/hr)
Control console(Operator Position)	
Outside patient entrance door	
Behind Windows (if applicable)	
Patient Waiting Area	

$$\text{Maximum Radiation level/week (mR/wk)} = \frac{\text{----- mAmin/week} \times \text{----Max. radiation level (mR/hr)}}{60 \times \text{-----mA used for measurement}}$$

Permissible limit

For location of Radiation Worker: 20 mSv in a year (40 mR/week)

For Location of Member of Public: 1 mSv in a year (2mR/week)